

BIO BM6+

Portable Digital Microscope Dermatoscope



www.portable-ultrasound-scanner.com

Quality, Professional, innovational

Portable Digital Microscope Dermatoscope BM6+


Features:

1. 3 inch TFT color display, rotated 360°
2. 8 LEDs with brightness adjustable
3. 5M sensor & professional microscopic lens
4. A metal stand fixes the dermatoscope
5. 20~230x magnification
6. Video, photo and sound capture
7. Works with Li-ion rechargeable battery or power supply
8. Micro SD card optional, download via USB
9. Measurement by software (1/1000mm)

Specifications:

Image sensor	5M pixels
Lens	High quality microscopic lens
Display	3 inches TFT panel,4:3 ratio
Magnification	20~230x
Still Image resolution	12M (4030x3024) ,9M,5M,3M,1.3M,VGA
Video capture resolution	VGA,QVGA
Frame rate	30fps@VGA
Light source	8 LEDs with brightness adjustable
Power source	(1) Li-ion 3.7V/1050mAh rechargeable battery (2) Adaptor:Input:100-240V;Output:5V 1A (3) USB
TV output	To any monitor with TV-in
Storage	Support Micro-SD card up to 32GB SDHC (not included)
OSD Language	English/Germany/French/Spanish/Italian/ Portuguese/Japanese/Chinese
Measurement	By Micro Capture Pro Software (Window XP/Vista/7/8,Mac 10.7 and above)
Software Language	English/Germany/French/Spanish
Dimension	130x103x29mm
Weight	195g(With battery)

Images:



Wuxi Biomedical Technology Co., Ltd. www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobile phone: 86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.